

Microliter 800 Series Syringes



The 800 series syringes are the reinforced plunger version of the original 700 series syringes. Attached to the flange of the syringe is a blue Syringe Holder. The Syringe Holder provides a spot to hold the syringe that will not heat up the barrel causing dispense inaccuracies. Additionally it provides support to small volume plungers that are prone to bending during injection. This series of syringes is ideal for dispensing volumes from 0.5 μ L up to 250 μ L. A variety of different terminations and needle options are available in this series.

Downloads	Technical Information	Contact Us +1-800-648-5950 sales@hamiltoncompany.com (mailto:sales@hamiltoncompany.com)
-----------	-----------------------	--

Displaying 1 - 18 of 18 products

Narrow Your Results	Part Number/REF	Description	Unit Price	Quantity
Volume	84852 (http://www.hamiltoncompany.com/products/syringes-and-needles/general-syringes/microliter-syringes/800-series/10-microL-Model-801-N-SYR-Cemented-NDL-26s-ga-2-in-point-style-2)	10 μL, Model 801 N SYR, Cemented NDL, 26s ga, 2 in, point style 2 (http://www.hamiltoncompany.com/products/syringes-and-needles/general-syringes/microliter-syringes/800-series/10-microL-Model-801-N-SYR-Cemented-NDL-26s-ga-2-in-point-style-2)	Quote Required	1 Add to Cart +
Termination	7642-01 (http://www.hamiltoncompany.com/products/syringes-and-needles/general-syringes/microliter-syringes/800-series/10-microL-Model-801-RN-SYR-NDL-Sold-Separately)	10 μL, Model 801 RN SYR, NDL Sold Separately (http://www.hamiltoncompany.com/products/syringes-and-needles/general-syringes/microliter-syringes/800-series/10-microL-Model-801-RN-SYR-NDL-Sold-Separately)	Quote Required	1 Add to Cart +
Gauge	84853 (http://www.hamiltoncompany.com/products/syringes-and-needles/general-syringes/microliter-syringes/800-series/10-microL-Model-801-RN-SYR-Small-Removable-NDL-26s-ga-2-in-point-style-2)	10 μL, Model 801 RN SYR, Small Removable NDL, 26s ga, 2 in, point style 2 (http://www.hamiltoncompany.com/products/syringes-and-needles/general-syringes/microliter-syringes/800-series/10-microL-Model-801-RN-SYR-Small-Removable-NDL-26s-ga-2-in-point-style-2)	Quote Required	1 Add to Cart +
Needle Length	84890 (http://www.hamiltoncompany.com/products/syringes-and-needles/general-syringes/microliter-syringes/800-series/10-25-50-100-and-250-microL-Model-RN-SYR-5-Pack-Removable-NDL-26s-and-22s-ga-2-in-point-style-2)	10, 25, 50, 100, and 250 μL, Model RN SYR, 5 Pack, Removable NDL, 26s and 22s ga, 2 in, point style 2 (http://www.hamiltoncompany.com/products/syringes-and-needles/general-syringes/microliter-syringes/800-series/10-25-50-100-and-250-microL-Model-RN-SYR-5-Pack-Removable-NDL-26s-and-22s-ga-2-in-point-style-2)	Quote Required	1 Add to Cart +
Point Style				

<p>84858 (http://www.hamiltoncompany.com/products/syringes-and-needles/general-syringes/microliter-syringes/800-series/100-microL-Model-810-N-SYR-Cemented-NDL-22s-ga-2-in-point-style-2)</p>	<p>100 µL, Model 810 N SYR, Cemented NDL, 22s ga, 2 in, point style 2 (http://www.hamiltoncompany.com/products/syringes-and-needles/general-syringes/microliter-syringes/800-series/100-microL-Model-810-N-SYR-Cemented-NDL-22s-ga-2-in-point-style-2)</p>	<p>Quote Required <input type="text" value="1"/></p>	<p>Add to Cart +</p>
<p>7645-01 (http://www.hamiltoncompany.com/products/syringes-and-needles/general-syringes/microliter-syringes/800-series/100-microL-Model-810-RN-SYR-NDL-Sold-Separately)</p>	<p>100 µL, Model 810 RN SYR, NDL Sold Separately (http://www.hamiltoncompany.com/products/syringes-and-needles/general-syringes/microliter-syringes/800-series/100-microL-Model-810-RN-SYR-NDL-Sold-Separately)</p>	<p>Quote Required <input type="text" value="1"/></p>	<p>Add to Cart +</p>
<p>84859 (http://www.hamiltoncompany.com/products/syringes-and-needles/general-syringes/microliter-syringes/800-series/100-microL-Model-810-RN-SYR-Small-Removable-NDL-22s-ga-2-in-point-style-2)</p>	<p>100 µL, Model 810 RN SYR, Small Removable NDL, 22s ga, 2 in, point style 2 (http://www.hamiltoncompany.com/products/syringes-and-needles/general-syringes/microliter-syringes/800-series/100-microL-Model-810-RN-SYR-Small-Removable-NDL-22s-ga-2-in-point-style-2)</p>	<p>Quote Required <input type="text" value="1"/></p>	<p>Add to Cart +</p>
<p>84854 (http://www.hamiltoncompany.com/products/syringes-and-needles/general-syringes/microliter-syringes/800-series/25-microL-Model-802-N-SYR-Cemented-NDL-22s-ga-2-in-point-style-2)</p>	<p>25 µL, Model 802 N SYR, Cemented NDL, 22s ga, 2 in, point style 2 (http://www.hamiltoncompany.com/products/syringes-and-needles/general-syringes/microliter-syringes/800-series/25-microL-Model-802-N-SYR-Cemented-NDL-22s-ga-2-in-point-style-2)</p>	<p>Quote Required <input type="text" value="1"/></p>	<p>Add to Cart +</p>
<p>7643-01 (http://www.hamiltoncompany.com/products/syringes-and-needles/general-syringes/microliter-syringes/800-series/25-microL-Model-802-RN-SYR-NDL-Sold-Separately)</p>	<p>25 µL, Model 802 RN SYR, NDL Sold Separately (http://www.hamiltoncompany.com/products/syringes-and-needles/general-syringes/microliter-syringes/800-series/25-microL-Model-802-RN-SYR-NDL-Sold-Separately)</p>	<p>Quote Required <input type="text" value="1"/></p>	<p>Add to Cart +</p>
<p>84855 (http://www.hamiltoncompany.com/products/syringes-and-needles/general-syringes/microliter-syringes/800-series/25-microL-Model-802-RN-SYR-Small-Removable-NDL-22s-ga-2-in-point-style-2)</p>	<p>25 µL, Model 802 RN SYR, Small Removable NDL, 22s ga, 2 in, point style 2 (http://www.hamiltoncompany.com/products/syringes-and-needles/general-syringes/microliter-syringes/800-series/25-microL-Model-802-RN-SYR-Small-Removable-NDL-22s-ga-2-in-point-style-2)</p>	<p>Quote Required <input type="text" value="1"/></p>	<p>Add to Cart +</p>
<p>84861 (http://www.hamiltoncompany.com/products/syringes-and-needles/general-syringes/microliter-syringes/800-series/250-microL-Model-825-RN-SYR-Large-Removable-NDL-22s-ga-2-in-point-style-2)</p>	<p>250 µL, Model 825 RN SYR, Large Removable NDL, 22s ga, 2 in, point style 2 (http://www.hamiltoncompany.com/products/syringes-and-needles/general-syringes/microliter-syringes/800-series/250-microL-Model-825-RN-SYR-Large-Removable-NDL-22s-ga-2-in-point-style-2)</p>	<p>Quote Required <input type="text" value="1"/></p>	<p>Add to Cart +</p>

<p>7646-01 (http://www.hamiltoncompany.com/products/syringes-and-needles/general-syringes/microliter-syringes/800-series/250-microL-Model-825-RN-SYR-NDL-Sold-Separately)</p>	<p>250 µL, Model 825 RN SYR, NDL Sold Separately (http://www.hamiltoncompany.com/products/syringes-and-needles/general-syringes/microliter-syringes/800-series/250-microL-Model-825-RN-SYR-NDL-Sold-Separately)</p>	<p>Quote Required <input type="text" value="1"/></p>	<p>Add to Cart +</p>
<p>84850 (http://www.hamiltoncompany.com/products/syringes-and-needles/general-syringes/microliter-syringes/800-series/5-microL-Model-85-N-SYR-Cemented-NDL-26s-ga-2-in-point-style-2)</p>	<p>5 µL, Model 85 N SYR, Cemented NDL, 26s ga, 2 in, point style 2 (http://www.hamiltoncompany.com/products/syringes-and-needles/general-syringes/microliter-syringes/800-series/5-microL-Model-85-N-SYR-Cemented-NDL-26s-ga-2-in-point-style-2)</p>	<p>Quote Required <input type="text" value="1"/></p>	<p>Add to Cart +</p>
<p>7641-01 (http://www.hamiltoncompany.com/products/syringes-and-needles/general-syringes/microliter-syringes/800-series/5-microL-Model-85-RN-SYR-NDL-Sold-Separately)</p>	<p>5 µL, Model 85 RN SYR, NDL Sold Separately (http://www.hamiltoncompany.com/products/syringes-and-needles/general-syringes/microliter-syringes/800-series/5-microL-Model-85-RN-SYR-NDL-Sold-Separately)</p>	<p>Quote Required <input type="text" value="1"/></p>	<p>Add to Cart +</p>
<p>84851 (http://www.hamiltoncompany.com/products/syringes-and-needles/general-syringes/microliter-syringes/800-series/5-microL-Model-85-RN-SYR-Small-Removable-NDL-26s-ga-2-in-point-style-2)</p>	<p>5 µL, Model 85 RN SYR, Small Removable NDL, 26s ga, 2 in, point style 2 (http://www.hamiltoncompany.com/products/syringes-and-needles/general-syringes/microliter-syringes/800-series/5-microL-Model-85-RN-SYR-Small-Removable-NDL-26s-ga-2-in-point-style-2)</p>	<p>Quote Required <input type="text" value="1"/></p>	<p>Add to Cart +</p>
<p>84856 (http://www.hamiltoncompany.com/products/syringes-and-needles/general-syringes/microliter-syringes/800-series/50-microL-Model-805-N-SYR-Cemented-NDL-22s-ga-2-in-point-style-2)</p>	<p>50 µL, Model 805 N SYR, Cemented NDL, 22s ga, 2 in, point style 2 (http://www.hamiltoncompany.com/products/syringes-and-needles/general-syringes/microliter-syringes/800-series/50-microL-Model-805-N-SYR-Cemented-NDL-22s-ga-2-in-point-style-2)</p>	<p>Quote Required <input type="text" value="1"/></p>	<p>Add to Cart +</p>
<p>7644-01 (http://www.hamiltoncompany.com/products/syringes-and-needles/general-syringes/microliter-syringes/800-series/50-microL-Model-805-RN-SYR-NDL-Sold-Separately)</p>	<p>50 µL, Model 805 RN SYR, NDL Sold Separately (http://www.hamiltoncompany.com/products/syringes-and-needles/general-syringes/microliter-syringes/800-series/50-microL-Model-805-RN-SYR-NDL-Sold-Separately)</p>	<p>Quote Required <input type="text" value="1"/></p>	<p>Add to Cart +</p>
<p>84857 (http://www.hamiltoncompany.com/products/syringes-and-needles/general-syringes/microliter-syringes/800-series/50-microL-Model-805-RN-SYR-Small-Removable-NDL-22s-ga-2-in-point-style-2)</p>	<p>50 µL, Model 805 RN SYR, Small Removable NDL, 22s ga, 2 in, point style 2 (http://www.hamiltoncompany.com/products/syringes-and-needles/general-syringes/microliter-syringes/800-series/50-microL-Model-805-RN-SYR-Small-Removable-NDL-22s-ga-2-in-point-style-2)</p>	<p>Quote Required <input type="text" value="1"/></p>	<p>Add to Cart +</p>

Display
[20](#) [50](#) [All](#)